

CH076 rabbit pAb

Cat No.:ES19874

For research use only

Overview

Product Name	CH076 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CH076 AA range: 149-199
Specificity	This antibody detects endogenous levels of human CH076
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Uncharacterized protein C8orf76
Gene Name	C8orf76
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84933
Human Swiss-Prot Number	Q96K31
Alternative Names	
Background	

