



# MRP4 rabbit pAb

Cat No.:ES19895

For research use only

## Overview

<b>Product Name</b>	MRP4 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB
<b>Species Cross-Reactivity</b>	Human;Rat;Mouse;
<b>Recommended dilutions</b>	WB 1:1000-2000
<b>Immunogen</b>	Synthesized peptide derived from human MRP4 AA range: 884-934
<b>Specificity</b>	This antibody detects endogenous levels of human MRP4
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C . Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Multidrug resistance-associated protein 4 (ATP-binding cassette sub-family C member 4) (MRP/cMOAT-related ABC transporter) (Multi-specific organic anion transporter B) (MOAT-B)
<b>Gene Name</b>	ABCC4 MRP4
<b>Cellular localization</b>	Basolateral cell membrane ; Multi-pass membrane protein . Apical cell membrane ; Multi-pass membrane protein . Its localization to the basolateral or apical membranes is tissue-dependent. .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	10257
<b>Human Swiss-Prot Number</b>	O15439
<b>Alternative Names</b>	
<b>Background</b>	

