



MFAP2 rabbit pAb

Cat No.:ES19903

For research use only

Overview

Product Name	MFAP2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MFAP2 AA range: 101-151
Specificity	This antibody detects endogenous levels of human MFAP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Microfibrillar-associated protein 2 (MFAP-2) (Microfibril-associated glycoprotein 1) (MAGP) (MAGP-1)
Gene Name	MFAP2 MAGP1
Cellular localization	Secreted, extracellular space, extracellular matrix.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	4237
Human Swiss-Prot Number	P55001
Alternative Names	
Background	

