



ISOC2 rabbit pAb

Cat No.:ES19904

For research use only

Overview

Product Name	ISOC2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ISOC2 AA range: 47-97
Specificity	This antibody detects endogenous levels of human ISOC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Isochorismatase domain-containing protein 2, mitochondrial
Gene Name	ISOC2
Cellular localization	Cytoplasm . Nucleus . Localizes to the nucleus in the presence of CDKN2A.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79763
Human Swiss-Prot Number	Q96AB3
Alternative Names	
Background	

