



EGFEM rabbit pAb

Cat No.:ES19914

For research use only

Overview

Product Name	EGFEM rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EGFEM AA range: 33-83
Specificity	This antibody detects endogenous levels of human EGFEM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Putative EGF-like and EMI domain-containing protein 1
Gene Name	EGFEM1P C3orf50 NCRNA00259
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q0D2K5
Alternative Names	
Background	

