



PHKG2 rabbit pAb

Cat No.:ES19924

For research use only

Overview

Product Name	PHKG2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PHKG2 AA range: 166-216
Specificity	This antibody detects endogenous levels of human PHKG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Phosphorylase b kinase gamma catalytic chain, testis/liver isoform (PHK-gamma-T) (EC 2.7.11.19) (PSK-C3) (Phosphorylase kinase subunit gamma-2) (Serine/threonine-protein kinase PHKG2) (EC 2.7.11.1) (E
Gene Name	PHKG2
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	5261
Human Swiss-Prot Number	P15735
Alternative Names	
Background	

