

GRD2I rabbit pAb

Cat No.:ES19933

For research use only

Overview

Product Name	GRD2I rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human GRD2I AA range: 43-93
Specificity	This antibody detects endogenous levels of human GRD2I
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Delphinin (Glutamate receptor, ionotropic, delta 2-interacting protein 1)
Gene Name	GRID2IP
Cellular localization	Cell junction, synapse, postsynaptic cell membrane .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	392862
Human Swiss-Prot Number	A4D2P6
Alternative Names	
Background	

