

ATG2B rabbit pAb

Cat No.:ES19935

For research use only

Overview

Product Name	ATG2B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ATG2B AA range: 266-316
Specificity	This antibody detects endogenous levels of human ATG2B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Autophagy-related protein 2 homolog B
Gene Name	ATG2B C14orf103
Cellular localization	Preautophagosomal structure membrane ; Peripheral membrane protein . Lipid droplet . Endoplasmic reticulum membrane ; Peripheral membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55102
Human Swiss-Prot Number	Q96BY7
Alternative Names	
Background	

