

DCNP1 rabbit pAb

Cat No.:ES19936

For research use only

Overview

Product Name	DCNP1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DCNP1 AA range: 27-77
Specificity	This antibody detects endogenous levels of human DCNP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Dendritic cell nuclear protein 1
Gene Name	DCNP1 C5orf20
Cellular localization	Nucleus . Cytoplasm . Particularly on the periphery. Colocalizes with corticotropin-releasing hormone (CRH) in parvocellular neurons in the paraventricular nucleus. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	140947
Human Swiss-Prot Number	Q8TF63
Alternative Names	
Background	

