



PHF2 rabbit pAb

Cat No.:ES19939

For research use only

Overview

| | |
|---------------------------------|---|
| Product Name | PHF2 rabbit pAb |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human;Mouse |
| Recommended dilutions | WB 1:1000-2000 |
| Immunogen | Synthesized peptide derived from human PHF2 AA range: 74-124 |
| Specificity | This antibody detects endogenous levels of human PHF2 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C . Avoid repeated freeze-thaw cycles. |
| Protein Name | Lysine-specific demethylase PHF2 (EC 1.14.11.-) (GRC5) (PHD finger protein 2) |
| Gene Name | PHF2 KIAA0662 |
| Cellular localization | Nucleus, nucleolus . Chromosome, centromere, kinetochore . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | 5253 |
| Human Swiss-Prot Number | O75151 |
| Alternative Names | |
| Background | |

