

CLD15 rabbit pAb

Cat No.:ES19940

For research use only

Overview

Product Name	CLD15 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CLD15 AA range: 66-116
Specificity	This antibody detects endogenous levels of human CLD15
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Claudin-15
Gene Name	CLDN15
Cellular localization	Cell junction, tight junction. Cell membrane; Multi-pass membrane protein. Tight junctions form continuous circumferential cell-cell contacts at the borders of apical and lateral cell membranes that seal the intercellular space and show up as strand-like structures in electron microscopy.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	24146
Human Swiss-Prot Number	P56746
Alternative Names	
Background	

