

CK049 rabbit pAb

Cat No.:ES19941

For research use only

Overview

Product Name	CK049 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CK049 AA range: 165-215
Specificity	This antibody detects endogenous levels of human CK049
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	UPF0705 protein C11orf49
Gene Name	C11orf49
Cellular localization	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriolar satellite . Cytoplasm, cytoskeleton . Associated with microtubules. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79096
Human Swiss-Prot Number	Q9H6J7
Alternative Names	
Background	

