



ENTP3 rabbit pAb

Cat No.:ES19944

For research use only

Overview

Product Name	ENTP3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ENTP3 AA range: 298-348
Specificity	This antibody detects endogenous levels of human ENTP3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Ectonucleoside triphosphate diphosphohydrolase 3 (NTPDase 3) (EC 3.6.1.5) (CD39 antigen-like 3) (Ecto-ATP diphosphohydrolase 3) (Ecto-ATPDase 3) (Ecto-ATPase 3) (Ecto-apyrase 3) (HB6)
Gene Name	ENTPD3 CD39L3
Cellular localization	Cell membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	956
Human Swiss-Prot Number	O75355
Alternative Names	
Background	

