



ITIH4 (80k, Cleaved-Arg689) rabbit pAb

Cat No.:ES19946

For research use only


Overview

Product Name	ITIH4 (80k, Cleaved-Arg689) rabbit pAb
Host species	Rabbit
Applications	WB;IHC
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300
Immunogen	Synthesized peptide derived from human ITIH4 (80k, Cleaved-Arg689)
Specificity	This antibody detects endogenous levels of Human ITIH4 (80k, Cleaved-Arg689)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	ITIH4 (80k, Cleaved-Arg689)
Gene Name	ITIH4 IHRP ITIHL1 PK120 PRO1851
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	80 100kD
Human Gene ID	3700
Human Swiss-Prot Number	Q14624
Alternative Names	Inter-alpha-trypsin inhibitor heavy chain H4 (ITI heavy chain H4;ITI-HC4;Inter-alpha-inhibitor heavy chain 4;Inter-alpha-trypsin inhibitor family heavy chain-related protein;IHRP;Plasma kallikrein sensitive glycoprotein 120;Gp120;PK-120) [Cleaved into: 70
Background	#REF!





ELK Biotechnology



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C