



# CATW (Cleaved-Val128) rabbit pAb

Cat No.:ES19963

For research use only

## Overview

<b>Product Name</b>	CATW (Cleaved-Val128) rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB; ELISA
<b>Species Cross-Reactivity</b>	Human;Mouse;Rat
<b>Recommended dilutions</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Immunogen</b>	Synthesized peptide derived from human CATW (Cleaved-Val128)
<b>Specificity</b>	This antibody detects endogenous levels of Human,Mouse,Rat CATW (Cleaved-Val128, protein was cleaved amino acid sequence between 127-128 )
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C . Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	CATW (Cleaved-Val128)
<b>Gene Name</b>	CTSW
<b>Cellular localization</b>	Endoplasmic reticulum .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	28 42kD
<b>Human Gene ID</b>	1521
<b>Human Swiss-Prot Number</b>	P56202
<b>Alternative Names</b>	Cathepsin W (EC 3.4.22.-;Lymphopain)
<b>Background</b>	

