

MYH-pan (Acetyl Lys1394) rabbit pAb

Cat No.:ES20118

For research use only

Overview

Product Name	MYH-pan (Acetyl Lys1394) rabbit pAb
Host species	Rabbit
Applications	WB;IHC
Species	Human;Mouse;Rat
Cross-Reactivity	
Recommended	WB 1:500-2000;IHC-p 1:50-300
dilutions	
Immunogen	Synthesized peptide derived from human MYH-pan (Acetyl Lys1394)
Specificity	This antibody detects endogenous levels of
	Human,Mouse,Rat MYH-pan (Acetyl Lys1394)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	MYH-pan (Acetyl Lys1394)
Gene Name	MYH1
Cellular localization	Cytoplasm, myofibril. Thick filaments of the myofibrils.
Purification	The antibody was affinity-purified from rabbit antiserum by
	affinity-chromatography using epitope-specific
	immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	77kD
Human Gene ID	4619
Human Swiss-Prot	P12882/Q9UKX2/P11055/Q9Y623/P13533/P12883/P13535
Number	
Alternative Names	Myosin-1 (Myosin heavy chain 1;Myosin heavy chain
	2x;MyHC-2x;Myosin heavy chain IIx/d;MyHC-IIx/d;Myosin
	heavy chain, skeletal muscle, adult 1)
Background	Myosin is a major contractile protein which converts
	chemical energy into mechanical energy through the
	hydrolysis of ATP. Myosin is a hexameric protein composed
	of a pair of myosin heavy chains (MYH) and two pairs of



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



nonidentical light chains. Myosin heavy chains are encoded by a multigene family. In mammals at least 10 different myosin heavy chain (MYH) isoforms have been described from striated, smooth, and nonmuscle cells. These isoforms show expression that is spatially and temporally regulated during development. [provided by RefSeq, Jul 2008],

Immunohistochemical analysis of paraffin-embedded human cervical carcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C