



NAGK (phospho-Ser74) rabbit pAb

Cat No.:ES20207

For research use only

Overview

Product Name	NAGK (phospho-Ser74) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human NAGK (phospho-Ser74)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat NAGK (phospho-Ser74)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	NAGK (phospho-Ser74)
Gene Name	NAGK
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	55577
Human Swiss-Prot Number	Q9UJ70
Alternative Names	N-acetyl-D-glucosamine kinase (N-acetylglucosamine kinase;EC 2.7.1.59;GlcNAc kinase)

Background

