

KITH_HHV1C rabbit pAb

Cat No.:ES20231

For research use only

Overview

Product Name	KITH_HHV1C rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from KITH_HHV1C
Specificity	This antibody detects endogenous levels of KITH_HHV1C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	KITH_HHV1C
Gene Name	KITH
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	P06478
Alternative Names	
Background	

