

# IL-17 rabbit pAb

Cat No.:ES20279

For research use only

## Overview

Product Name	IL-17 rabbit pAb
Host species	Rabbit
Applications	WB;Flow Cyt; ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000;Flow Cyt 1:50-200; ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human IL-17 AA range: 50-100
Specificity	This antibody detects endogenous levels of Human IL-17
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	IL-17
Gene Name	IL17A CTLA8 IL17
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3605
Human Swiss-Prot Number	Q16552
Alternative Names	Interleukin-17A (IL-17;IL-17A;Cytotoxic T-lymphocyte-associated antigen 8;CTLA-8)
Background	function:Induces stromal cells to produce proinflammatory and hematopoietic cytokines. Enhances the surface expression of the intracellular adhesion molecule-1 (ICAM-1) in fibroblasts.,online information:Interleukin-17 entry,PTM:Found both in glycosylated and nonglycosylated forms.,similarity:Belongs to the IL-17





**ELK Biotechnology**

family,subunit:Homodimer,tissue  
specificity:Restricted to activated memory T-cells.,



+86-27-59760950

ELKbio@ELKbiotech.com

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C