



TSLP rabbit pAb

Cat No.:ES20289

For research use only

Overview

Product Name	TSLP rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human TSLP AA range: 91-140
Specificity	This antibody detects endogenous levels of Human TSLP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	TSLP
Gene Name	TSLP
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	85480
Human Swiss-Prot Number	Q969D9
Alternative Names	Thymic stromal lymphopoietin
Background	function:Cytokine that induces the release of T cell-attracting chemokines from monocytes and, in particular, enhances the maturation of CD11c(+) dendritic cells. Can induce allergic inflammation by directly activating mast cells.,induction:Released by primary epithelial cells in response to certain microbial products, physical injury, or inflammatory cytokines.,subunit:Interacts with a receptor composed of CRLF2 and IL7R.,tissue





ELK Biotechnology

specificity: Expressed in a number of tissues including heart, liver and prostate.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C