



Calponin 1 rabbit pAb

Cat No.:ES20296

For research use only

Overview

Product Name	Calponin 1 rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human Calponin 1 AA range: 111-160
Specificity	This antibody detects endogenous levels of Human Calponin 1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Calponin 1
Gene Name	CNN1
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	1264
Human Swiss-Prot Number	P51911
Alternative Names	Calponin-1 (Basic calponin;Calponin H1, smooth muscle)
Background	

