

## AGRP rabbit pAb

## Cat No.:ES20302

For research use only

## Overview

Product Name	AGRP rabbit pAb	
Host species	Rabbit	
Applications	WB; ELISA	
Species Cross-Reactivity	Human;Rat;Mouse;	
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000	
Immunogen	Synthesized peptide derived from human AGRP AA	
	range: 21-70	
Specificity	This antibody detects endogenous levels of Human	
	AGRP	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20 $^\circ\!\mathrm{C}$ . Avoid repeated freeze-thaw cycles.	
Protein Name	AGRP	
Gene Name	AGRP AGRT ART	
Cellular localization	Secreted . Golgi apparatus lumen .	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band		
Human Gene ID	181	
Human Swiss-Prot Number	O00253	
Alternative Names	Agouti-related protein	
Background	disease:Defects in AGRP may be a cause of	
	autosomal dominant obesity	
	[MIM:601665].,domain:The presence of a 'disulfide	
	through disulfide knot' structurally defines this	
	protein as a knottin.,function:Plays a role in weight	
	homeostasis. May play a role in the regulation of	al
	melanocortin receptors within the hypothalamus	
	and adrenal gland, and therefore in the central	
	control of feeding., similarity: Contains 1 agouti	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



domain.,tissue specificity:Expressed primarily in the adrenal gland, subthalamic nucleus, and hypothalamus, with a lower level of expression occurring in testis, lung, and kidney.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C