

UL33_HHV6U rabbit pAb

Cat No.:ES20314

For research use only

Overview

Product Name	UL33_HHV6U rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human virus
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human UL33_HHV6U
Specificity	This antibody detects endogenous levels of Human herpesvirus 6A UL33_HHV6U
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	UL33_HHV6U
Gene Name	UL33_HHV6U
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	P52380
Alternative Names	
Background	

