



IGF1 rabbit pAb

Cat No.:ES20322

For research use only

Overview

Product Name	IGF1 rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human IGF1
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat IGF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	IGF1
Gene Name	IGF1 IBP1
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3479
Human Swiss-Prot Number	P05019
Alternative Names	Insulin-like growth factor I (IGF-I;Mechano growth factor;MGF;Somatomedin-C)

Background

