



GRAB rabbit pAb

Cat No.:ES20325

For research use only

Overview

Product Name	GRAB rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Mouse
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from mouse GRAB
Specificity	This antibody detects endogenous levels of Mouse GRAB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	GRAB
Gene Name	GRAB
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	P04187
Alternative Names	Granzyme B (EC 3.4.21.79;C11;CTLA-1;Cathepsin G-like 1;CTSL1;Cytotoxic T-lymphocyte proteinase 2;Lymphocyte protease;Fragmentin-2;Granzyme-2;Human lymphocyte protein;HLP;SECT;T-cell serine protease 1-3E)

Background

