



IL17F rabbit pAb

Cat No.:ES20327

For research use only

Overview

Product Name	IL17F rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Mouse
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from mouse IL17F
Specificity	This antibody detects endogenous levels of Mouse IL17F
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	IL17F
Gene Name	IL17F
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q7TNI7
Alternative Names	Interleukin-17F (IL-17F;Cytokine ML-1;Interleukin-24;IL-24)
Background	

