



IL27A rabbit pAb

Cat No.:ES20329

For research use only

Overview

| | |
|---------------------------------|---|
| Product Name | IL27A rabbit pAb |
| Host species | Rabbit |
| Applications | WB; ELISA |
| Species Cross-Reactivity | Mouse |
| Recommended dilutions | WB 1:1000-2000 ELISA 1:5000-20000 |
| Immunogen | Synthesized peptide derived from mouse IL27A |
| Specificity | This antibody detects endogenous levels of Mouse IL27A |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Protein Name | IL27A |
| Gene Name | IL27A |
| Cellular localization | |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | |
| Human Gene ID | |
| Human Swiss-Prot Number | Q8K3I6 |
| Alternative Names | Interleukin-27 subunit alpha (IL-27 subunit alpha;IL-27-A;IL27-A;p28) |
| Background | |

