

# IL5 rabbit pAb

Cat No.:ES20331

For research use only

## Overview

Product Name	IL5 rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from N-terminal mouse IL5
Specificity	This antibody detects endogenous levels of Mouse,Rat IL5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	IL5
Gene Name	IL5
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q08125/P04401
Alternative Names	
Background	

