



FGF7 rabbit pAb

Cat No.:ES20352

For research use only

Overview

Product Name	FGF7 rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human FGF7
Specificity	This antibody detects endogenous levels of Human,Rat FGF7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	FGF7
Gene Name	FGF7 KGF
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	2252
Human Swiss-Prot Number	P21781
Alternative Names	Fibroblast growth factor 7 (FGF-7;Heparin-binding growth factor 7;HBGF-7;Keratinocyte growth factor)
Background	

