



PHGDH rabbit pAb

Cat No.:ES20402

For research use only

Overview

Product Name	PHGDH rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human PHGDH
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat PHGDH
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	PHGDH
Gene Name	PHGDH PGDH3
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	26227
Human Swiss-Prot Number	O43175
Alternative Names	D-3-phosphoglycerate dehydrogenase (3-PGDH;EC 1.1.1.95)
Background	

