

IL-9 rabbit pAb

Cat No.:ES20413

For research use only

Overview

Product Name	IL-9 rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human IL-9
Specificity	This antibody detects endogenous levels of Human IL-9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	IL-9
Gene Name	IL9
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3578
Human Swiss-Prot Number	P15248
Alternative Names	Interleukin-9 (IL-9;Cytokine P40;T-cell growth factor P40)
Background	function:Supports IL-2 independent and IL-4 independent growth of helper T-cells.,online information:Interleukin-9 entry,similarity:Belongs to the IL-7/IL-9 family.,

