

IL-3 rabbit pAb

Cat No.:ES20425

For research use only

Overview

Product Name	IL-3 rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human IL-3
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat IL-3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	IL-3
Gene Name	IL3
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3562
Human Swiss-Prot Number	P08700
Alternative Names	Interleukin-3 (IL-3;Hematopoietic growth factor;Mast cell growth factor;MCGF;Multipotential colony-stimulating factor;P-cell-stimulating factor)

Background

