

IL-7 rabbit pAb

Cat No.:ES20434

For research use only

Overview

Product Name	IL-7 rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human IL-7
Specificity	This antibody detects endogenous levels of Human IL-7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	IL-7
Gene Name	IL7
Cellular localization	Secreted.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3574
Human Swiss-Prot Number	P13232
Alternative Names	Interleukin-7 (IL-7)
Background	function:Hematopoietic growth factor capable of stimulating the proliferation of lymphoid progenitors. It is important for proliferation during certain stages of B-cell maturation.,online information:Interleukin-7 entry,similarity:Belongs to the IL-7/IL-9 family.,

