

ZBTB45 rabbit pAb


Cat No.:ES20595

For research use only

Overview

Product Name	ZBTB45 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB: 1:500-1000 IHC: 1:200-500
Immunogen	Recombinant Protein of ZBTB45
Specificity	The antibody detects endogenous ZBTB45 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger and BTB domain-containing protein 45
Gene Name	ZBTB45
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	
Observed band	70kD
Human Gene ID	84878
Human Swiss-Prot Number	Q96K62
Alternative Names	Zinc finger and BTB domain-containing protein 45; Zinc finger protein 499
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 4 C2H2-type zinc fingers.,





Western blot analysis of 1) Hela, 2) Jurkat, 3) Mouse Brain, 4) Mouse Kidney, 5) Rat Brain using ZBTB45 Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000

Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using ZBTB45 Polyclonal Antibody.

