



GABA A Receptor γ 2 rabbit pAb

Cat No.:ES20747

For research use only

Overview


Product Name	GABA A Receptor γ 2 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF
Species Cross-Reactivity	Human;Rat;Mouse
Recommended dilutions	WB 1:1000-2000, IHC 1:100-200
Immunogen	Synthetic Peptide of GABA A Receptor γ 2
Specificity	GABA A Receptor γ 2 protein(A229) detects endogenous levels of GABA A Receptor γ 2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Gamma-aminobutyric acid receptor subunit gamma-2 (GABA(A) receptor subunit gamma-2)
Gene Name	GABRG2
Cellular localization	Cell junction, synapse, postsynaptic cell membrane ; Multi-pass membrane protein . Cell membrane ; Multi-pass membrane protein . Cell projection, dendrite . Cytoplasmic vesicle membrane .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kD
Human Gene ID	2566
Human Swiss-Prot Number	P18507
Alternative Names	Gamma-aminobutyric acid receptor subunit gamma-2 (GABA(A) receptor subunit gamma-2)
Background	This gene encodes a gamma-aminobutyric acid (GABA) receptor. GABA is the major inhibitory neurotransmitter in the mammalian brain, where it acts at GABA-A receptors, which are ligand-gated chloride channels. GABA-A receptors are






pentameric, consisting of proteins from several subunit classes: alpha, beta, gamma, delta and rho. Mutations in this gene have been associated with epilepsy and febrile seizures. Multiple transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jul 2008],

Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue with GABA A Receptor $\gamma 2$ Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using GABA A Receptor $\gamma 2$ Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using GABA A Receptor $\gamma 2$ Rabbit pAb diluted at 1:200.

