



GABA A Receptor α 2 rabbit pAb

Cat No.:ES20750

For research use only

Overview


| | |
|---------------------------------|--|
| Product Name | GABA A Receptor α 2 rabbit pAb |
| Host species | Rabbit |
| Applications | WB;IHC;IF |
| Species Cross-Reactivity | Human;Rat;Mouse |
| Recommended dilutions | WB 1:1000-2000, IHC 1:100-200 |
| Immunogen | Synthetic Peptide of GABA A Receptor α 2 |
| Specificity | GABA A Receptor α 2 protein(A224) detects endogenous levels of GABA A Receptor α 2 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Protein Name | Gamma-aminobutyric acid receptor subunit alpha-2 (GABA(A) receptor subunit alpha-2) |
| Gene Name | GABRA2 |
| Cellular localization | Cell junction, synapse, postsynaptic cell membrane ; Multi-pass membrane protein . Cell membrane ; Multi-pass membrane protein . Cytoplasmic vesicle membrane . Cell projection, dendrite . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 51kD |
| Human Gene ID | 282236 |
| Human Swiss-Prot Number | P47869 |
| Alternative Names | Gamma-aminobutyric acid receptor subunit alpha-2 (GABA(A) receptor subunit alpha-2) |
| Background | GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA-A receptors, which are ligand-gated chloride channels. Chloride conductance of these channels can be modulated by agents such as benzodiazepines that bind to the |





GABA-A receptor. At least 16 distinct subunits of GABA-A receptors have been identified. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2013],

Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue with GABA A Receptor $\alpha 2$ Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using GABA A Receptor $\alpha 2$ Rabbit pAb diluted at 1:200.

