



eEF1A2 binding protein rabbit pAb

Cat No.:ES20768

For research use only

Overview


Product Name	eEF1A2 binding protein rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF
Species Cross-Reactivity	Human;Rat;Mouse;Pig
Recommended dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200
Immunogen	Recombinant Protein of eEF1A2 binding protein
Specificity	The antibody detects endogenous eEF1A2 binding protein protein
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Immunoglobulin-like and fibronectin type III domain-containing protein 1 (EEF1A2-binding protein 1) (KY-interacting protein 1)
Gene Name	IGFN1
Cellular localization	Nucleus . Cytoplasm, myofibril, sarcomere, Z line .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	115kD
Human Gene ID	91156
Human Swiss-Prot Number	Q86VF2
Alternative Names	Immunoglobulin-like and fibronectin type III domain-containing protein 1 (EEF1A2-binding protein 1;KY-interacting protein 1)
Background	similarity:Contains 4 fibronectin type-III domains.,similarity:Contains 5 Ig-like (immunoglobulin-like) domains.,subunit:Interacts with KY.,tissue specificity:Expressed in skeletal muscle.,





ELK Biotechnology

Western blot analysis of Pig Skeletal Muscle with eEF1A2
Rabbit pAb diluted at 1:1,000



Immunohistochemical analysis of paraffin-embedded Pig
Skeletal Muscle using eEF1A2Rabbit pAb diluted at 1:100



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C