



ECFP rabbit pAb

Cat No.:ES20769

For research use only

Overview

Product Name	ECFP rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Species independent
Recommended dilutions	WB: 1:2000-5000
Immunogen	Recombinant Protein of ECFP
Specificity	The antibody detects ECFP & CFP tag protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	
Gene Name	
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	
Observed band	
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	YM3228
Background	Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish <i>Aequorea victoria</i> .





ELK Biotechnology



Western blot analysis of ECFP recombinant protein,
diluted at 1:5000. Secondary antibody(catalog#:RS0002)
was diluted at 1:20000



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C