



ELK Biotechnology

# Actin beta rabbit pAb(Caenorhabditis elegans)

Cat No.:ES20846

For research use only

---

## Overview

|                          |   |
|--------------------------|---|
| Product Name             | Actin beta rabbit pAb(Caenorhabditis elegans)   |
| Host species             | Rabbit  |
| Applications             | WB  |
| Species Cross-Reactivity | C.elegans   |
| Recommended dilutions    | WB 1:5000-20000   |
| Immunogen                | Purified recombinant Full length Actin beta (Caenorhabditis elegans)  |
| Specificity              | This antibody detects endogenous levels of Actin beta (Caenorhabditis elegans)  |
| Formulation              | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| Storage                  | Store at -20°C. Avoid repeated freeze-thaw cycles.  |
| Protein Name             | Actin beta  |
| Gene Name                |   |
| Cellular localization    |   |
| Purification             | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality                | Polyclonal  |
| Concentration            | 1 mg/ml   |
| Observed band            |   |
| Human Gene ID            |   |
| Human Swiss-Prot Number  |   |
| Alternative Names        | Beta Actin, $\beta$ Actin   |
| Background               |   |



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei, P.R.C.



**ELK Biotechnology**



Western Blot analysis of *C. elegans* cells using Actin beta rabbit pAb(*Caenorhabditis elegans*). Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

ELKbio@ELKbiotech.com

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C