

GAPDH rabbit pAb(Caenorhabditis elegans)

Cat No.:ES20847

For research use only


Overview

Product Name	GAPDH rabbit pAb(Caenorhabditis elegans)
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	C.elegans
Recommended dilutions	WB 1:5000-20000
Immunogen	Purified recombinant Full length GAPDH (Caenorhabditis elegans)
Specificity	This antibody detects endogenous levels of GAPDH (Caenorhabditis elegans)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	GAPDH
Gene Name	
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	P04970
Alternative Names	GAPDH
Background	





ELK Biotechnology



Western Blot analysis of *C. elegans* cells using GAPDH rabbit pAb(*Caenorhabditis elegans*). Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C