

# GAPDH rabbit pAb(Drosophilidae)

Cat No.:ES20852

For research use only


## Overview

Product Name	GAPDH rabbit pAb(Drosophilidae)
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Drosophilidae
Recommended dilutions	WB 1:5000-20000
Immunogen	Purified recombinant Full length GAPDH (Drosophilidae)
Specificity	This antibody detects endogenous levels of GAPDH (Drosophilidae)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	GAPDH
Gene Name	
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	P07486
Alternative Names	GAPDH
Background	





**ELK Biotechnology**



Western Blot analysis of Drosophilidae cells using GAPDH rabbit pAb(Drosophilidae). Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

ELKbio@ELKbiotech.com

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C