



# Tubulin alpha rabbit pAb(Yeast)

Cat No.:ES20854

For research use only


## Overview

<b>Product Name</b>	Tubulin alpha rabbit pAb(Yeast)
<b>Host species</b>	Rabbit
<b>Applications</b>	WB
<b>Species Cross-Reactivity</b>	Yeast
<b>Recommended dilutions</b>	WB 1:5000-20000
<b>Immunogen</b>	Purified recombinant Full length Tubulin alpha (Yeast)
<b>Specificity</b>	This antibody detects endogenous levels of Tubulin alpha (Yeast)
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Tubulin alpha
<b>Gene Name</b>	
<b>Cellular localization</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	
<b>Human Swiss-Prot Number</b>	
<b>Alternative Names</b>	Tubulin alpha
<b>Background</b>	





**ELK Biotechnology**



Western Blot analysis of Yeast cells using Tubulin alpha rabbit pAb(Yeast). Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

[ELKbio@ELKbiotech.com](mailto:ELKbio@ELKbiotech.com)

[www.elkbiotech.com](http://www.elkbiotech.com)

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C