

## CYP26A1 rabbit pAb

**Cat No.: ES2095** 

For research use only

## Overview

Product Name CYP26A1 rabbit pAb

Host species Rabbit
Applications WB;IF;ELISA

Species Cross-Reactivity Human; Mouse; Rat

**Recommended dilutions** Western Blot: 1/500 - 1/2000. Immunofluorescence:

1/200 - 1/1000. ELISA: 1/20000. Not yet tested in

other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human Cytochrome P450

26A1. AA range:251-300

**Specificity** CYP26A1 Polyclonal Antibody detects endogenous

levels of CYP26A1 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Cytochrome P450 26A1

Gene Name CYP26A1

**Cellular localization** Endoplasmic reticulum membrane ; Peripheral

membrane protein. Microsome membrane;

Peripheral membrane protein.

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 56kD
Human Gene ID 1592
Human Swiss-Prot Number O43174

Alternative Names CYP26A1; CYP26; P450RAI1; Cytochrome P450

26A1; Cytochrome P450 retinoic acid-inactivating 1; Cytochrome P450RAI; hP450RAI; Retinoic acid 4-hydroxylase; Retinoic acid-metabolizing

cytochrome

+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com






**Background** 

cytochrome P450 family 26 subfamily A member 1(CYP26A1) Homo sapiens This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This endoplasmic reticulum protein acts on retinoids, including all-trans-retinoic acid (RA), with both 4-hydroxylation and 18-hydroxylation activities. This enzyme regulates the cellular level of retinoic acid which is involved in regulation of gene expression in both embryonic and adult tissues. Two alternatively spliced transcript variants of this gene, which encode the distinct isoforms, have been reported. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells using CYP26A1 Polyclonal Antibody diluted at 1:1000




Western Blot analysis of 293 cells using CYP26A1 Polyclonal Antibody diluted at 1:1000



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





Western blot analysis of lysates from HepG2, HeLa, and 293 cells, using Cytochrome P450 26A1 Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950