

DAPK3 rabbit pAb

Cat No.:ES2143

For research use only

Overview

Product Name	DAPK3 rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/20000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
-	peptide derived from human DAPK3. AA
	range:241-290
Specificity	DAPK3 Polyclonal Antibody detects endogenous
	levels of DAPK3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Death-associated protein kinase 3
Gene Name	DAPK3
Cellular localization	Nucleus . Cytoplasm . Predominantly localizes to the
	cytoplasm but can shuttle between the nucleus and
	cytoplasm; cytoplasmic localization is promoted by
	phosphorylation at Thr-299 and involves Rho/Rock
	signaling; [Isoform 1]: Nucleus . Cytoplasm .; [I
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	52kD
Human Gene ID	1613
Human Swiss-Prot Number	O43293
Alternative Names	DAPK3; ZIPK; Death-associated protein kinase 3;
	DAP kinase 3; DAP-like kinase; Dlk; MYPT1 kinase;
	Zipper-interacting protein kinase; ZIP-kinase



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

(kD)

117-85-

48-

34-


26-


19-

HuvEc

Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells using DAPK3 Polyclonal Antibody





Western blot analysis of lysates from HUVEC cells, using DAPK3 Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C