

DNA pol λ rabbit pAb

Cat No.: ES2186

For research use only

Overview

Product Name DNA pol λ rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications. The antiserum was produced against synthesized

Immunogen The antiserum was produced against synthesized peptide derived from human DNA Polymerase

lambda. AA range:451-500

Specificity DNA pol λ Polyclonal Antibody detects endogenous

levels of DNA pol λ protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name DNA polymerase lambda

Gene Name POLL

Cellular localization Nucleus . Chromosome . Accumulates at sites of

DNA damage. .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 63kD
Human Gene ID 27343
Human Swiss-Prot Number Q9UGP5

Alternative Names POLL; DNA polymerase lambda; Pol Lambda; DNA

polymerase beta-2; Pol beta2; DNA polymerase

kappa

Background This gene encodes a DNA polymerase. DNA


polymerases catalyze DNA-template-directed extension of the 3'-end of a DNA strand. This




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com




particular polymerase, which is a member of the X family of DNA polymerases, likely plays a role in non-homologous end joining and other DNA repair processes. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2010],



Western Blot analysis of various cells using DNA pol λ Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absor



Western blot analysis of lysates from 293 cells, using DNA Polymerase lambda Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950 EL

ELKbio@ELKbiotech.com

www.elkbiotech.com