

FASTKD1 rabbit pAb

Cat No.:ES2330

For research use only

Overview

Product Name	FASTKD1 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human; Mouse; Monkey
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not
	yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human FAKD1. AA
	range:561-610
Specificity	FASTKD1 Polyclonal Antibody detects endogenous
	levels of FASTKD1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	FAST kinase domain-containing protein 1
Gene Name	FASTKD1
Cellular localization	Mitochondrion . Preferentially localizes to
	mitochondrial RNA granules, platforms for
	post-transcriptional RNA modification and ribosome
	assembly (PubMed:28335001)
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	97kD
Human Gene ID	79675
Human Swiss-Prot Number	Q53R41
Alternative Names	FASTKD1; KIAA1800; FAST kinase domain-containing
	protein 1
Background	caution:It is uncertain whether Met-1 or Met-15 is
	the initiator., similarity: Belongs to the FAST kinase
	family., similarity: Contains 1 RAP domain.,




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C







ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C