

Granzyme K rabbit pAb

Cat No.: ES2470

For research use only

Overview

Product Name Granzyme K rabbit pAb

Host species Rabbit

Applications WB;IHC;IF;ELISA Species Cross-Reactivity Human;Mouse;Rat

Recommended dilutions Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA:

1/20000. Not yet tested in other applications.
The antiserum was produced against synthesized

Immunogen The antiserum was produced against synthesized peptide derived from human GRAK. AA

range:61-110

Specificity Granzyme K Polyclonal Antibody detects

endogenous levels of Granzyme K protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Granzyme K

Gene Name GZMK

Cellular localization Secreted. Cytoplasmic granule.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 33kD
Human Gene ID 3003
Human Swiss-Prot Number P49863

Alternative Names GZMK; TRYP2; Granzyme K; Fragmentin-3;

Granzyme-3; NK-tryptase-2; NK-Tryp-2

Background This gene product is a member of a group of related

serine proteases from the cytoplasmic granules of cytotoxic lymphocytes. Cytolytic T lymphocytes (CTL) and natural killer (NK) cells share the




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.




remarkable ability to recognize, bind, and lyse specific target cells. They are thought to protect their host by lysing cells bearing on their surface 'nonself' antigens, usually peptides or proteins resulting from infection by intracellular pathogens. The protein described here lacks consensus sequences for N-glycosylation present in other granzymes. [provided by RefSeq, Jul 2008],


Western Blot analysis of various cells using Granzyme K Polyclonal Antibody diluted at 1:1000



Immunofluorescence analysis of NIH/3T3 cells, using GRAK Antibody. The picture on the right is blocked with the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using GRAK Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950







Western blot analysis of lysates from Jurkat cells, using GRAK Antibody. The lane on the right is blocked with the synthesized peptide.

