

## **HGK** rabbit pAb

Cat No.: ES2511

For research use only

## Overview

Product Name HGK rabbit pAb

Host species Rabbit

Applications WB;IF;ELISA Species Cross-Reactivity Human;Mouse

**Recommended dilutions** Western Blot: 1/500 - 1/2000. Immunofluorescence:

1/200 - 1/1000. ELISA: 1/40000. Not yet tested in

other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human MEKKK 4. AA

range:406-455

**Specificity** HGK Polyclonal Antibody detects endogenous levels

of HGK protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Mitogen-activated protein kinase kinase

kinase 4

**Gene Name** MAP4K4 **Cellular localization** Cytoplasm .

**Purification** The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 142kD
Human Gene ID 9448
Human Swiss-Prot Number 095819

Alternative Names MAP4K4; HGK; KIAA0687; NIK; Mitogen-activated

protein kinase kinase kinase kinase 4; HPK/GCK-like kinase HGK; MAPK/ERK kinase kinase kinase 4; MEK kinase kinase 4; MEKKK 4; Nck-interacting kinase mitogen-activated protein kinase kinase kinase

kinase 4(MAP4K4) Homo sapiens The protein




Background

+86-27-59760950 ELKbio@ELKbiotech.com ww



encoded by this gene is a member of the serine/threonine protein kinase family. This kinase has been shown to specifically activate MAPK8/JNK. The activation of MAPK8 by this kinase is found to be inhibited by the dominant-negative mutants of MAP3K7/TAK1, MAP2K4/MKK4, and MAP2K7/MKK7, which suggests that this kinase may function through the MAP3K7-MAP2K4-MAP2K7 kinase cascade, and mediate the TNF-alpha signaling pathway. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],


170-130-95-72-55Western Blot analysis of various cells using HGK Polyclonal Antibody



Immunofluorescence analysis of A549 cells, using MEKKK 4 Antibody. The picture on the right is blocked with the synthesized peptide.







Western blot analysis of lysates from HUVEC and COLO cells, using MEKKK 4 Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950