

LPD lipase rabbit pAb

Cat No.:ES2723

For research use only

Overview

Product Name	LPD lipase rabbit pAb
Host species	Rabbit
Applications	WB;IHC;IF;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	Western Blot: 1/500 - 1/2000.
	Immunohistochemistry: 1/100 - 1/300. ELISA:
	1/10000. Not yet tested in other applications.
Immunogen	The antiserum was produced against synthesized
-	peptide derived from human LIPI. AA range:289-338
Specificity	LPD lipase Polyclonal Antibody detects endogenous
	levels of LPD lipase protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Lipase member I
Gene Name	LIPI
Cellular localization	[Isoform 1]: Cell membrane . Secreted . May
	associate with lipid draft; [Isoform 2]: Cell
	membrane . Secreted . May associate with lipid
	draft
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	53kD
Human Gene ID	149998
Human Swiss-Prot Number	Q6XZB0
Alternative Names	LIPI; LPDL; PRED5; Lipase member I; LIPI;
	Cancer/testis antigen 17; CT17; LPD lipase;
	Membrane-associated phosphatidic acid-selective
	phospholipase A1-beta; mPA-PLA1 beta
Background	The protein encoded by this gene is a phospholipase



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



COLO205

(kD)

117-85-

48-

34-

26-

19-

LIPI--

(kD)

117-

85-

48-

34-

26-

19-

COLO COLO

-- 117

-- 85

-- 48

-- 34

-- 26 -- 19 (kD)

LIPI

that hydrolyzes phosphatidic acid to produce lysophosphatidic acid. Defects in this gene are a cause of susceptibility to familial hypertrigliceridemia. This gene is also expressed at high levels in Ewing family tumor cells. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2014],

Western Blot analysis of various cells using LPD lipase Polyclonal Antibody

Western blot analysis of lysates from COLO cells, using LIPI Antibody. The lane on the right is blocked with the synthesized peptide.

Western blot analysis of the lysates from HepG2 cells using LIPI antibody.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C