

MKP-1 rabbit pAb

Cat No.: ES2795

For research use only

Overview

Product Name MKP-1 rabbit pAb

Host species Rabbit

WB;IHC;IF;ELISA **Applications Species Cross-Reactivity** Human; Rat; Mouse;

Recommended dilutions Western Blot: 1/500 - 1/2000.

> Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.

The antiserum was produced against synthesized **Immunogen**

peptide derived from human MKP1. AA

range:318-367

Specificity MKP-1 Polyclonal Antibody detects endogenous

levels of MKP-1 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Store at -20°C. Avoid repeated freeze-thaw cycles. **Storage**

Dual specificity protein phosphatase 1 **Protein Name**

Gene Name DUSP1 Cellular localization Nucleus.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml **Observed band** 39kD **Human Gene ID** 1843 **Human Swiss-Prot Number** P28562

Alternative Names DUSP1; CL100; MKP1; PTPN10; VH1; Dual specificity

> protein phosphatase 1; Dual specificity protein phosphatase hVH1; Mitogen-activated protein kinase phosphatase 1; MAP kinase phosphatase 1;

MKP-1; Protein-tyrosine phosphatase CL100

The expression of DUSP1 gene is induced in human **Background**


skin fibroblasts by oxidative/heat stress and growth




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com



factors. It specifies a protein with structural features similar to members of the non-receptor-type protein-tyrosine phosphatase family, and which has significant amino-acid sequence similarity to a Tyr/Ser-protein phosphatase encoded by the late gene H1 of vaccinia virus. The bacterially expressed and purified DUSP1 protein has intrinsic phosphatase activity, and specifically inactivates mitogen-activated protein (MAP) kinase in vitro by the concomitant dephosphorylation of both its phosphothreonine and phosphotyrosine residues. Furthermore, it suppresses the activation of MAP kinase by oncogenic ras in extracts of Xenopus oocytes. Thus, DUSP1 may play an important role in the human cellular response to environmental stress as well as in the negative regulation of cellular proliferati



Western Blot analysis of various cells using MKP-1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).



Western Blot analysis of Jurkat cells using MKP-1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003,Inventbiotech,MN,USA).




+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com





Immunohistochemistry analysis of paraffin-embedded human brain tissue, using MKP1 Antibody. The picture on the right is blocked with the synthesized peptide.



+86-27-59760950

Western blot analysis of lysates from Jurkat cells, using MKP1 Antibody. The lane on the right is blocked with the synthesized peptide.

